

	SK-5000A	SK-5000B
Measurement mode	Fully automatic/Automatic/Manual	
Corneal thickness measurement	Not support	Support Measurement range-CCT:400-800um Measurement increment-CCT:10um Light Source Blue LED
Display	7-inch 1024*600 color touch screen, support up and down 180°, left and right 60° flip	
Fixation light	1 internal fixation light (green LED), 2 external fixation lights	
Measurement accuracy	1mmHg (0.133Kpa)	
Measurement range	0-60mmHg (0-8Kpa) Include automatic switching between 30/60mmHg	
Printer	Built-in thermal printer, automatic paper cutting	
Data output/input	Input: USB; Output:LAN/RS232C	
Size and weight	500*300*500, <25kg	



Chongqing Sunkingdom Medical Instrument Co.,Ltd.

Consulting Tel: : 023-68102787
After sales support: : 17316790887
Official website: www.cqsunkingdom.com



FULLY AUTOMATIC NON-CONTACT TONOMETER

QUICK SPEED
3D AUTO-TRACKING
AUTO MEASUREMENT
SOFT AIR PUFF, SAFETY,
PACHYETER+TONOMETER



Sun Kingdom 上邦医疗
Chongqing Sunkingdom Medical Instrument Co.,Ltd.

Non-contact tonometer
SK-5000A/B



Electric chin rest

SK-5000B supports corneal thickness measurement
More functions




Auto alignment
Auto focus

Fully automatic mode
One touch
You can get OS&OD data
It also supports semi-automatic/manual mode





Corneal thickness measurement
Auto soft air puff

Support air puff detection after turn on
Support safety limit setting
Safely equipment use and operation



Auto switch OS&OD



Print report
Auto paper cutting



ONE CLICK

FULLY
AUTOMATIC



7 inch color capacitive touch screen
Can be flip up and down 180 °
Flip right and left 60 °
Operation and observation are more convenient

